



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MTG16 Antibody [#C10912]

Catalog Number: C10912

Concentration: 1mg/ml

Swiss-Prot No.: O75081/Q0P5Z6

Other Names: Myeloid translocation gene on chromosome 16 protein; hMTG16; MTG8-related protein 2; Zinc finger MYND domain-containing protein 4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MTG16.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MTG16 antibody detects endogenous levels of total MTG16 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Gamou T., Blood 91:4028-4037(1998).

Calabi F., Genomics 52:332-341(1998).

Salomon-Nguyen F., Leukemia 14:1704-1705(2000).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using MTG16 antibody.