



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TEF Antibody [#C10916]

Catalog Number: C10916

Concentration: 1mg/ml

Swiss-Prot No.: Q10587

Other Names: thyrotroph embryonic factor; thyrotrophic embryonic factor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TEF.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TEF antibody detects endogenous levels of total TEF protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Khatib Z.A., Genomics 23:344-351(1994).

Hunger S.P., Blood 87:4607-4617(1996).

Collins J.E., Genome Biol. 5:RESEARCH84.1-RESEARCH84.11(2004).

