



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AKAP13 Antibody [#C10921]

Catalog Number: C10921

Concentration: 1mg/ml

Swiss-Prot No.: Q12802

Other Names: AKAP 13, Protein kinase A-anchoring protein 13, Breast cancer nuclear receptor-binding auxiliary protein, Human thyroid-anchoring protein 31, Guanine nucleotide exchange factor Lbc, AKAP-Lbc, LBC oncogene

P47, Lymphoid blast crisis oncogene, Non-oncogenic Rho GTPase-specific GTP exchange factor, BRX, HT31, LBC

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human AKAP13.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AKAP13 antibody detects endogenous levels of total AKAP13 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:20000

References:

Toksoz D., Oncogene 9:621-628(1994).


Rubino D., Oncogene 16:2513-2526(1998).

Sterpetti P., Mol. Cell. Biol. 19:1334-1345(1999).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using AKAP13 antibody.



Western blot analysis of extracts from Jurkat cells and A549 cells, using AKAP13 antibody.