



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

Nuclear Factor 1 Antibody [#C10925]

Catalog Number: C10925

Concentration: 1mg/ml

Swiss-Prot No.: Q12857/O00712/P08651/Q14938

Other Names: CCAAT-box binding transcription factor, CTF, KIAA1439, NF-I/A, NF1-A, NFIA, Nuclear factor 1 A-type, Nuclear factor 1/A, TGGCA-binding protein

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

Nuclear Factor 1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Nuclear Factor 1 antibody detects endogenous levels of total Nuclear Factor 1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Nagase T., DNA Res. 7:65-73(2000).

Gregory S.G., Nature 441:315-321(2006).

Qian F., Genomics 28:66-73(1995).

