



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ELAVL2 Antibody [#C10931]

Catalog Number: C10931

Concentration: 1mg/ml

Swiss-Prot No.: Q12926

Other Names: Hu-antigen B, HuB, ELAV-like neuronal protein 1, Nervous system-specific RNA-binding protein Hel-N1, HUB

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ELAVL2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ELAVL2 antibody detects endogenous levels of total ELAVL2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

King P.H., J. Neurosci. 14:1943-1952(1994).

Gao F.B., Proc. Natl. Acad. Sci. U.S.A. 91:11207-11211(1994).

Humphray S.J., Nature 429:369-374(2004).

