



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FOXC1/2 Antibody [#C10937]

Catalog Number: C10937

Concentration: 1mg/ml

Swiss-Prot No.: Q12948

Other Names: Forkhead-related protein, FKHL7,
Forkhead-related transcription factor 3, FREAC-3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FOXC1/2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FOXC1/2 antibody detects endogenous levels of total FOXC1/2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Nishimura D.Y., Nat. Genet. 19:140-147(1998).

Mears A.J., Am. J. Hum. Genet. 63:1316-1328(1998).

Mungall A.J., Nature 425:805-811(2003).

