



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FOXL1 Antibody [#C10939]

Catalog Number: C10939

Concentration: 1mg/ml

Swiss-Prot No.: Q12952

Other Names: ILF, ILF-1, ILF1, cellular transcription factor ILF-1, forkhead box K2, interleukin enhancer-binding 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FOXL1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FOXL1 antibody detects endogenous levels of total FOXL1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Fang J., Am. J. Hum. Genet. 67:1382-1388(2000).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Pierrou S., EMBO J. 13:5002-5012(1994).

