



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CAF1A Antibody [#C10945]

Catalog Number: C10945

Concentration: 1mg/ml

Swiss-Prot No.: Q13111

Other Names: Chromatin assembly factor 1 subunit A, CAF-1 subunit A, Chromatin assembly factor I p150 subunit, CAF-I 150 kDa subunit, CAF-I p150, hp150

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CAF1A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: CAF1A antibody detects endogenous levels of total CAF1A protein.

Reactivity: Human; Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Dong H., Gene 264:187-196(2001) [PubMed: 11250073].
Grimwood J., Nature 428:529-535(2004) [PubMed: 15057824].
Kaufman P.D., Cell 81:1105-1114(1995) [PubMed: 7600578].



For Research Use Only