



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CAF1B Antibody [#C10946]

Catalog Number: C10946

Concentration: 1mg/ml

Swiss-Prot No.: Q13112

Other Names: CAF-1 subunit B, Chromatin assembly factor I p60 subunit, CAF-I 60 kDa subunit, CAF-Ip60, M-phase phosphoprotein 7

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

CAF1B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CAF1B antibody detects endogenous levels of total CAF1B protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Kaufman P.D., Cell 81:1105-1114(1995).

Hattori M., Nature 405:311-319(2000).

Matsumoto-Taniura N., Mol. Biol. Cell 7:1455-1469(1996).

