



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

CGK 2 Antibody [#C10961]

Catalog Number: C10961

Concentration: 1mg/ml

Swiss-Prot No.: Q13237

Other Names: cGMP-dependent protein kinase 2, EC 2.7.11.12, CGK 2, Cgkii, Type II cGMP-dependent protein kinase, PRKG2, PRKGR2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CGK 2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CGK 2 antibody detects endogenous levels of total CGK 2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Orstavik S., Biochem. Biophys. Res. Commun. 220:759-765(1996).

Fujii M., FEBS Lett. 375:263-267(1995).

Witczak O., Biochem. Biophys. Res. Commun. 245:113-119(1998).

