



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## SFRS5 Antibody [#C10962]

**Catalog Number:** C10962

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q13243

**Other Names:** Splicing factor, arginine/serine-rich 5, Pre-mRNA-splicing factor SRP40, Delayed-early protein HRS, HRS, SRP40

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human SFRS5.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** SFRS5 antibody detects endogenous levels of total SFRS5 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 100%, Positives = 100%);  
Rat (Identities = 100%, Positives = 100%)

**Applications:** IF: 1:500~1:1000 ELISA: 1:20000

**References:**

Screaton G.R., EMBO J. 14:4336-4349(1995).  
Du K., Gene 204:243-249(1997).  
Zahler A.M., Genes Dev. 6:837-847(1992).

