



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## MTA1 Antibody [#C10969]

**Catalog Number:** C10969

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q13330

**Other Names:** Metastasis-associated protein MTA1, MTA1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human MTA1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MTA1 antibody detects endogenous levels of total MTA1 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Toh Y., J. Biol. Chem. 269:22958-22963(1994).

Kumar R., Nature 418:654-657(2002).

Xue Y., Mol. Cell 2:851-861(1998).

