



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KLF Antibody [#C10971]

Catalog Number: C10971

Concentration: 1mg/ml

Swiss-Prot No.: Q13351

Other Names: Krueppel-like factor 1, Erythroid krueppel-like transcription factor, EKLF, Erythroid transcription factor, EKLF

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human KLF.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KLF antibody detects endogenous levels of total KLF protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Bieker J.J., DNA Cell Biol. 15:347-352(1996).

van Ree J.H., Genomics 39:393-395(1997).

Grimwood J., Nature 428:529-535(2004).

