



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FOXE3 Antibody [#C10979]

Catalog Number: C10979

Concentration: 1mg/ml

Swiss-Prot No.: Q13461

Other Names: Forkhead box protein E3,
Forkhead-related protein FKHL12, Forkhead-related
transcription factor 8, FREAC-8, FOXE3, FKHL12,
FREAC8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

FOXE3.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: FOXE3 antibody detects endogenous levels
of total FOXE3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Semina E.V., Hum. Mol. Genet. 10:231-236(2001).

Gregory S.G., Nature 441:315-321(2006).

Larsson C., Genomics 30:464-469(1995).

