



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Smad4 Antibody [#C10983]

Catalog Number: C10983

Concentration: 1mg/ml

Swiss-Prot No.: Q13485

Other Names: SMAD 4, Mothers against decapentaplegic homolog 4, Mothers against DPP homolog 4, Deletion target in pancreatic carcinoma 4, hSMAD4, DPC4, MADH4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

Smad4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Smad4 antibody detects endogenous levels of total Smad4 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Hahn S.A., Science 271:350-353(1996).

Zhang Y., Nature 383:168-172(1996).

Moskaluk C.A., Diagn. Mol. Pathol. 6:85-90(1997).



For Research Use Only