



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TDG Antibody [#C10996]

Catalog Number: C10996

Concentration: 1mg/ml

Swiss-Prot No.: Q13569

Other Names: C-JUN leucine zipper interactive protein JZA-3; EC 3.2.2.-; G/T mismatch-specific thymine DNA glycosylase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TDG.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TDG antibody detects endogenous levels of total TDG protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Neddermann P., J. Biol. Chem. 271:12767-12774(1996).
The MGC Project Team, Genome Res. 14:2121-2127(2004).
Hardeland U., EMBO J. 21:1456-1464(2002).

