



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DGKZ Antibody [#C10999]

Catalog Number: C10999

Concentration: 1mg/ml

Swiss-Prot No.: Q13574

Other Names: EC 2.7.1.107; Diglyceride kinase zeta; DGK-zeta; DAG kinase zeta

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human DGKZ.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DGKZ antibody detects endogenous levels of total DGKZ protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Bunting M., J. Biol. Chem. 271:10230-10236(1996).

Ding L., Proc. Natl. Acad. Sci. U.S.A. 94:5519-5524(1997).

Topham M.K., Nature 394:697-700(1998).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using DGKZ antibody.