



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TF2H2 Antibody [#C11010]

Catalog Number: C11010

Concentration: 1mg/ml

Swiss-Prot No.: Q13888

Other Names: TFIIH basal transcription factor complex p44 subunit, Basic transcription factor 2 44 kDa subunit, BTF2-p44, General transcription factor IIH polypeptide 2, GTF2H2, BTF2P44

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human

TF2H2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TF2H2 antibody detects endogenous levels of total TF2H2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Humbert S., EMBO J. 13:2393-2398(1994).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

van der Steege G., Eur. J. Hum. Genet. 3:87-95(1995).

