



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

GIT2 Antibody [#C11019]

Catalog Number: C11019

Concentration: 1mg/ml

Swiss-Prot No.: Q14161

Other Names: ARF GTPase-activating protein GIT2; G protein-coupled receptor kinase- interactor 2; KIAA0148; Tyrosine-phosphorylated protein CAT-2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GIT2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GIT2 antibody detects endogenous levels of total GIT2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Premont R.T., J. Biol. Chem. 275:22373-22380(2000).

Nagase T., DNA Res. 2:167-174(1995).

Nakajima D., DNA Res. 9:99-106(2002).

