



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

E2F2 Antibody [#C11025]

Catalog Number: C11025

Concentration: 1mg/ml

Swiss-Prot No.: Q14209

Other Names: Transcription factor E2F2, E2F-2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human E2F2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: E2F2 antibody detects endogenous levels of total E2F2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Ivey-Hoyle M., Mol. Cell. Biol. 13:7802-7812(1993).

Gregory S.G., Nature 441:315-321(2006).

Novy M., Mol. Biol. Cell 16:2181-2190(2005).

