



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

Glucokinase Regulator Antibody [#C11030]

Catalog Number: C11030

Concentration: 1mg/ml

Swiss-Prot No.: Q14397

Other Names: Glucokinase regulatory protein,
Glucokinase regulator

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Glucokinase Regulator.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Glucokinase Regulator antibody detects endogenous levels of total Glucokinase Regulator protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Warner J.P., Mamm. Genome 6:532-536(1995).

Hayward B.E., Genomics 49:137-142(1998).

Rieder M.J., Submitted (JUN-2003) to the EMBL/GenBank/DDBJ databases.

