



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## GLPK3 Antibody [#C11031]

**Catalog Number:** C11031

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q14409

**Other Names:** Glycerol kinase, testis specific 1, EC 2.7.1.30, ATP:glycerol 3-phosphotransferase,

Glycerokinase, GK, GK3P, GKP3, GKTB

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from internal of human GLPK3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** GLPK3 antibody detects endogenous levels of total GLPK3 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Sargent C.A., Hum. Mol. Genet. 3:1317-1324(1994).

