



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

GLPK2 Antibody [#C11032]

Catalog Number: C11032

Concentration: 1mg/ml

Swiss-Prot No.: Q14410

Other Names: Glycerol kinase, testis specific 2, EC 2.7.1.30, ATP:glycerol 3-phosphotransferase,

Glycerokinase, GK, GK2, GKP2, GKTA

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human GLPK2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GLPK2 antibody detects endogenous levels of total GLPK2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Sargent C.A., Hum. Mol. Genet. 3:1317-1324(1994).

