



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAD4 Antibody [#C11039]

Catalog Number: C11039

Concentration: 1mg/ml

Swiss-Prot No.: Q14582

Other Names: Max-associated protein 4; MAX dimerization protein 4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MAD4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAD4 antibody detects endogenous levels of total MAD4 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:10000

References:

Hillier L.W., Nature 434:724-731(2005).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Pribill I., Submitted (JAN-1998) to the EMBL/GenBank/DDBJ databases.



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using MAD4 antibody.