



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

NR1I3 Antibody [#C11053]

Catalog Number: C11053

Concentration: 1mg/ml

Swiss-Prot No.: Q14994

Other Names: CAR, Constitutive androstane receptor, Orphan nuclear receptor MB67, Orphan nuclear receptor NR1I3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

NR1I3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: NR1I3 antibody detects endogenous levels of total NR1I3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Baes M., Mol. Cell. Biol. 14:1544-1552(1994).

Choi H.S., J. Steroid Biochem. Mol. Biol. 56:23-30(1996).

Xu R.X., Mol. Cell 16:919-928(2004).

