



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF76 Antibody [#C10573]

Catalog Number: C10573

Concentration: 1mg/ml

Swiss-Prot No.: P36508

Other Names: Zinc finger protein 76, ZNF76, D6S229E

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ZNF76.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF76 antibody detects endogenous levels of total ZNF76 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:10000

References:

Ragoussis J., Genomics 14:673-679(1992).

Phelan M., Submitted (OCT-2004) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using ZNF76 antibody.