



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

hnRNP G Antibody [#C10581]

Catalog Number: C10581

Concentration: 1mg/ml

Swiss-Prot No.: P38159

Other Names: HNRPG, RBMX, RNA binding motif protein, X chromosome, ROG, glycoprotein p43, heterogeneous nuclear ribonucleoprotein G

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human hnRNP G.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: hnRNP G antibody detects endogenous levels of total hnRNP G protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Soulard M., Nucleic Acids Res. 21:4210-4217(1993).

Ross M.T., Nature 434:325-337(2005).

Venables J.P., Hum. Mol. Genet. 8:959-969(1999).



Western blot analysis of extracts from HeLa cells, using hnRNP G antibody.



Immunofluorescence analysis of HepG2 cells, using hnRNP G antibody.