



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ETV6 Antibody [#C10593]

Catalog Number: C10593

Concentration: 1mg/ml

Swiss-Prot No.: P41212

Other Names: ETS translocation variant 6, ETS-related protein Tel1, ETV6, TEL1, Transcription factor ETV6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ETV6.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ETV6 antibody detects endogenous levels of total ETV6 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:1000

References:


Golub T.R., Cell 77:307-316(1994).

Baens M., Genome Res. 6:404-413(1996).

Mullighan C.G., Nature 446:758-764(2007).



Western blot analysis of extracts from HeLa cells, HepG2 cells and Jurkat cells, using ETV6 antibody.



Immunofluorescence analysis of HeLa cells, using ETV6 antibody.