



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BUD31 Antibody [#C10594]

Catalog Number: C10594

Concentration: 1mg/ml

Swiss-Prot No.: P41223

Other Names: BUD31, EDG-2, G10

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human BUD31.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BUD31 antibody detects endogenous levels of total BUD31 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Hla T., Biochim. Biophys. Acta 1260:227-229(1995).

Hillier L.W., Nature 424:157-164(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

