



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

STAT5A Antibody [#C10605]

Catalog Number: C10605

Concentration: 1mg/ml

Swiss-Prot No.: P42229

Other Names: MGF, MPF, Mammary gland factor, STA5A, STAT5, STAT5A, Signal transducer and activator of transcription 5A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

STAT5A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: STAT5A antibody detects endogenous levels of total STAT5A protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Hou J., Immunity 2:321-329(1995).

Litterst C.M., J. Biol. Chem. 278:45340-45351(2003).

Martens N., J. Biol. Chem. 280:13817-13823(2005).

