



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## GATA4 Antibody [#C10618]

**Catalog Number:** C10618

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P43694

**Other Names:** DNA binding protein GATA-GT2, GAT4, GATA binding factor-4, GATA-4, Transcription factor GATA-4

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

GATA4.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** GATA4 antibody detects endogenous levels of total GATA4 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Huang W.Y., Gene 155:219-223(1995).

Garg V., Nature 424:443-447(2003).

