



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CBX5 Antibody [#C10620]

Catalog Number: C10620

Concentration: 1mg/ml

Swiss-Prot No.: P45973

Other Names: HP1 alpha, HP1A, antigen p25, chromobox protein homolog 5, heterochromatin protein 1 alpha

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

CBX5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CBX5 antibody detects endogenous levels of total CBX5 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Saunders W.S., J. Cell Sci. 104:573-582(1993).

Ye Q., J. Biol. Chem. 271:14653-14656(1996).

Minc E., Chromosoma 108:220-234(1999).

