



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ATRX Antibody [#C10627]

Catalog Number: C10627

Concentration: 1mg/ml

Swiss-Prot No.: P46100

Other Names: EC 3.6.1.-, RAD54L, Transcriptional regulator ATRX, X-linked helicase II, X-linked nuclear protein, XH2, XNP, Znf-HX

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human ATRX.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ATRXantibody detects endogenous levels of total ATRX protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Picketts D.J., Hum. Mol. Genet. 5:1899-1907(1996).

Villard L., Genomics 43:149-155(1997).

Kitano T., Mol. Biol. Evol. 20:1281-1289(2003).

