



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAP2K3 Antibody [#C10629]

Catalog Number: C10629

Concentration: 1mg/ml

Swiss-Prot No.: P46734

Other Names: Dual specificity mitogen-activated protein kinase kinase 3; MAP kinase kinase 3; MAP2K3; MAPK/ERK kinase 3; MAPKK 3; MAPKK3; MEK 3; MEK3; mitogen-activated protein kinase kinase 3; MKK3; MP2K3; PRKMK3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from N-terminal of human MAP2K3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAP2K3 antibody detects endogenous levels of total MAP2K3 protein.


Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Derijard B., Science 267:682-685(1995).



For Research Use Only