



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ID4 Antibody [#C10631]

**Catalog Number:** C10631

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P47928

**Other Names:** DNA-binding protein inhibitor ID-4,  
Inhibitor of DNA binding 4, ID4

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ID4.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ID4 antibody detects endogenous levels of total ID4 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Rigolet M., DNA Res. 5:309-313(1998).

Pagliuca A., Genomics 27:200-203(1995).

Mungall A.J., Nature 425:805-811(2003).

