



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TISD Antibody [#C10632]

Catalog Number: C10632

Concentration: 1mg/ml

Swiss-Prot No.: P47974

Other Names: Butyrate response factor 2, Protein TIS11D, EGF-response factor 2, ERF-2, ZFP36L2, BRF2, ERF2, RNF162C, TIS11D

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TISD.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TISD antibody detects endogenous levels of total TISD protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:13000 ELISA: 1:20000

References:

Ino T., *Oncogene* 11:2705-2710(1995).

Nie X.F., *Gene* 152:285-286(1995).

Hudson B.P., *Nat. Struct. Mol. Biol.* 11:257-264(2004).

