



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDK8 Antibody [#C10643]

Catalog Number: C10643

Concentration: 1mg/ml

Swiss-Prot No.: P49336

Other Names: Cell division protein kinase 8, EC 2.7.11.22, EC 2.7.11.23, Protein kinase K35, CDK8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human CDK8.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDK8 antibody detects endogenous levels of total CDK8 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Tassan J.-P., Proc. Natl. Acad. Sci. U.S.A. 92:8871-8875(1995).

Fryer C.J., Mol. Cell 16:509-520(2004).

