



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

CEBPE Antibody [#C10650]

Catalog Number: C10650

Concentration: 1mg/ml

Swiss-Prot No.: P49716

Other Names: C/EBP delta

Nuclear factor NF-IL6-beta

NF-IL6-beta

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CEBPE.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CEBPE antibody detects endogenous levels of total CEBPE protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Kinoshita S., Proc. Natl. Acad. Sci. U.S.A. 89:1473-1476(1992).

Cleutjens C.B., Genomics 16:520-523(1993).

Wang J.-M., Nucleic Acids Res. 34:217-231(2006).

