



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## p57KIP2 Antibody [#C10659]

**Catalog Number:** C10659

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P49918

**Other Names:** CDKN1C, CDN1C, CDNC,  
Cyclin-dependent kinase inhibitor 1C, Cyclin-dependent  
kinase inhibitor p57, KIP2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered  
saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl,  
0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against  
synthesized peptide derived from internal of human

p57KIP2.

**Purification:** The antibody was affinity-purified from  
rabbit antiserum by affinity-chromatography using  
epitope-specific immunogen.

**Specificity:** p57KIP2 antibody detects endogenous levels  
of total p57KIP2 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)


**Applications:** IHC: 1:50~1:100 ELISA: 1:10000

**References:**

Matsuoka S., Genes Dev. 9:650-662(1995).

Reid L.H., Cancer Res. 56:1214-1218(1996).

Tokino T., Hum. Genet. 97:625-631(1996).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded  
human brain tissue, using p57KIP2 antibody.

**For Research Use Only**