



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## MEOX2 Antibody [#C10662]

**Catalog Number:** C10662

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P50222

**Other Names:** Mesenchyme homeobox 2

Growth arrest-specific homeobox

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human MEOX2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MEOX2 antibody detects endogenous levels of total MEOX2 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Grigoriou M., Genomics 26:550-555(1995).

Lepage D.F., Genomics 24:535-540(1994).

Hillier L.W., Nature 424:157-164(2003).

