



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ER81 Antibody [#C10667]

Catalog Number: C10667

Concentration: 1mg/ml

Swiss-Prot No.: P50549

Other Names: ER81, ETS translocation variant 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human ER81.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ER81 antibody detects endogenous levels of total ER81 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Jeon I.-S., Oncogene 10:1229-1234(1995).

Monte D., Oncogene 11:771-779(1995).

Hillier L.W., Nature 424:157-164(2003).

