



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## RORG Antibody [#C10676]

**Catalog Number:** C10676

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P51449

**Other Names:** Nuclear receptor ROR-gamma;  
Retinoid-related orphan receptor-gamma; Nuclear  
receptor RZR-gamma; Nuclear receptor subfamily 1  
group F member 3; RORC; NR1F3; RORG; RZRG

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered  
saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl,  
0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against  
synthesized peptide derived from internal of human  
RORG.

**Purification:** The antibody was affinity-purified from  
rabbit antiserum by affinity-chromatography using  
epitope-specific immunogen.

**Specificity:** RORG antibody detects endogenous levels  
of total RORG protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 100%, Positives = 100%)

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Hirose T., Biochem. Biophys. Res. Commun.  
205:1976-1983(1994).  
Bechtel S., BMC Genomics 8:399-399(2007).  
Gregory S.G., Nature 441:315-321(2006).

