



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

GALK1 Antibody [#C10680]

Catalog Number: C10680

Concentration: 1mg/ml

Swiss-Prot No.: P51570

Other Names: Galactokinase, EC 2.7.1.6,
Galactose kinase, GALK

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human GALK1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GALK1 antibody detects endogenous levels of total GALK1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Stambolian D., Nat. Genet. 10:307-312(1995).

Bergsma D.J., Genome Res. 6:980-985(1996).

Ai Y., Biochem. Biophys. Res. Commun. 212:687-691(1995).

