



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BRCA2 Antibody [#C10682]

Catalog Number: C10682

Concentration: 1mg/ml

Swiss-Prot No.: P51587

Other Names: Breast cancer type 2 susceptibility protein

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human BRCA2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BRCA2 antibody detects endogenous levels of total BRCA2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Rat (Identities = 84%, Positives = 92%)


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Wooster R., Nature 378:789-792(1995).

Tavtigian S.V., Nat. Genet. 12:333-337(1996).

Dunham A., Nature 428:522-528(2004).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using BRCA2 antibody.