



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAP3KL4 Antibody [#C11190]

Catalog Number: C11190

Concentration: 1mg/ml

Swiss-Prot No.: Q5TCX8

Other Names: Mitogen-activated protein kinase kinase, EC 2.7.11.25, Mixed lineage kinase 4, MLK4, KIAA1804

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MAP3KL4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAP3KL4 antibody detects endogenous levels of total MAP3KL4 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Kvasha S., Submitted (APR-2001) to the EMBL/GenBank/DDBJ databases.

Nagase T., DNA Res. 8:85-95(2001).

Gregory S.G., Nature 441:315-321(2006).

