



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LDOCL Antibody [#C11219]

Catalog Number: C11219

Concentration: 1mg/ml

Swiss-Prot No.: Q6ICCC9

Other Names: Protein LDOC1L, Leucine zipper protein down-regulated in cancer cells-like, Mammalian retrotransposon-derived protein 6, LDOC1L, MAR6, MART6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

LDOCL.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LDOCL antibody detects endogenous levels of total LDOCL protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Collins J.E., Genome Biol. 5:RESEARCH84.1-RESEARCH84.11(2004).
Dunham I., Nature 402:489-495(1999).
Ota T., Nat. Genet. 36:40-45(2004).

