



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## SFXN4 Antibody [#C11227]

**Catalog Number:** C11227

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q6P4A7

**Other Names:** Breast cancer resistance marker 1,  
Sideroflexin-4, BCRM1, SFXN4

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human SFXN4.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** SFXN4 antibody detects endogenous levels of total SFXN4 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Zheng H., DNA Seq. 14:369-373(2003).

Deloukas P., Nature 429:375-381(2004).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

