



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## MOL2C Antibody [#C11242]

**Catalog Number:** C11242

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q70IA8

**Other Names:** Mps one binder kinase activator-like 2C, Mob1 homolog 2C, Protein Mob3C, MOBKL2C, MOB3C

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human MOL2C.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MOL2C antibody detects endogenous levels of total MOL2C protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Florindo C.S., Submitted (AUG-2003) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team, omc Res. 14:2121-2127(2004).

