



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

MAEA Antibody [#C11253]

Catalog Number: C11253

Concentration: 1mg/ml

Swiss-Prot No.: Q7L5Y9

Other Names: Macrophage erythroblast attacher, Erythroblast macrophage protein, Cell proliferation-inducing gene 5 protein, Lung cancer-related protein 10, HLC-10, MAEA, EMP, HLC10, PIG5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

MAEA.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAEA antibody detects endogenous levels of total MAEA protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Hanspal M., Blood 92:2940-2950(1998).

Bala S., Biochem. Biophys. Res. Commun. 342:1040-1048(2006).

