



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

STRAD Antibody [#C11258]

Catalog Number: C11258

Concentration: 1mg/ml

Swiss-Prot No.: Q7RTN6

Other Names: STE20-related adapter protein, STRAD alpha, Serologically defined breast cancer antigen NY-BR-96, STRAD, LYK5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

STRAD.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: STRAD antibody detects endogenous levels of total STRAD protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Scanlan M.J., Cancer Immun. 1:4-4(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Baas A.F., EMBO J. 22:3062-3072(2003).

